


PERFORM-X

Advanced High Frequency Radiographic Systems

Simplicity, Performance and Safety



The Perform-X radiographic systems are Control-X Medical's most advanced radiographic systems due to the various component configurations which can be made to match your diagnostic imaging needs.

A perfect choice for any radiographic environment, take advantage of these system's versatility, durability and hospital grade quality performance when searching for top quality radiographic systems.

ADVANCED HIGH FREQUENCY RADIOGRAPHIC SYSTEMS

GENERATORS	SPECIFICATIONS
Available models	Milestone HF LC Series, MS HF STR 100 Series, PCX HF Series
kW range (depending on model)	32 / 40 / 50 / 65 / 80
kVp range (depending on model)	40 - 150
mA range	300 / 400 / 500 / 650 / 800 / 1000

- Microprocessor controlled generators with over 100 APR programs
- Simple and user-friendly APR technique reprogramming procedures

TUBESTANDS	SPECIFICATIONS
Available models	TS 99, TS 99T, TS 05R, TS 05RT, Halo OTS
Movements (depending on model)	Vertical, longitudinal, rotational column (R), transverse arm (T)
X-ray tube rotation	± 180°
Mounting options (depending on model)	Floor, floor-to-wall, floor-to-ceiling, ceiling

- Heavy-duty, hospital grade quality performance
- Extended vertical track; electromagnetic brake system and dual rope counterbalancing system

TABLES	SPECIFICATIONS
Available models	Stylix 2, Phoenix
Movements	Four-way float top (Stylix 2 and Phoenix), elevating (Phoenix)
Weight limit	Useful use of 500 lbs (230 kg), Capable use of 660 lbs (300 kg)
Controls	All tabletop movements are controlled by conveniently located foot pedals

- Precision bearing system allows the tabletop to glide effortlessly
- Electromagnetic brake system for the tabletop
- DR compatible

WALLSTANDS	SPECIFICATIONS
Available models	WS 99, WS 99T
Movements (depending on model)	Vertical, tilting (T)
Mounting	Floor-to-wall, floor

- Durability: designed to withstand the rigors of the heaviest of caseloads
- Extended vertical track; electromagnetic brake system and dual rope counterbalancing system
- DR compatible

For more details on each component, please refer to their respective brochure.

For more information, please visit our website at www.cxmed.com or contact your local authorized distributor.



Control-X Medical
 1755 Atlas Street
 Columbus, OH 43228, USA
 Tel: 1-614-777-9729
 Toll free: 1-800-777-9729
 Fax: 1-614-777-0395
 E-mail: sales@cxmed.com
 Web: www.cxmed.com



Advanced Diagnostic Services
 4374 State Street, Bettendorf, IA 52722 | phone (563) 275-7211
 E-mail: scotton@adsqc.com

www.cxmed.com